

**KEEP YOUR CODE CARD IN A SAFE PLACE. NO CODE? NO RESULT!**



You can ask for the result by using the unique code on the code card.

**SCAN ME**



You can ask for the result at <https://testalize.me/uitslag> within 48 hours. To that end, use the code on the code card.

You can use this sampling kit to test for Chlamydia and Gonorrhoea. Your material will be examined at our laboratory (ISO15189). The results are

## RESULT

## STD-TEST

# INSTRUCTIONS

### STD-TEST ORAL

Chlamydia and Gonorrhoea



**MAIL THE RETURN BOX**

You can take the test as late as 14 days after you had (unsafe) sex for the last time.

# CONTENTS OF THE BOX



return box



instructions



code card



cotton swab\*



vial\*\*



sealbag\*\*\*

\* cotton swab: Aptima® Multitest Swab (SW) CE0845 MDD 93/42/EEC | Hologic Inc.

\*\* vial: Aptima® Multitest Swab Transport Medium (STM) | Hologic Inc.

\*\*\* sealbag Daklapack

## HOW TO CONDUCT THE TEST

Leave the fluid in the vial and avoid skin contact. Do not ingest. Testing is simple. There are only four steps to complete:

- 1 Wash your hands. Remove the cotton swab from its packaging. Do not touch the cotton and make sure it is not exposed to dirt.
- 2 Hold the cotton swab on the snap line and carefully place it in your mouth. Move it along your tonsils (both sides of your mouth) and the back of your throat. Take the cotton swab out of your mouth without touching your tongue or cheeks.
- 3 Open the vial. Place the cotton swab inside and snap it along the line. Securely close the vial. Place the vial in the sealbag and seal the bag. Place the full sealbag in the return box and close the box.
- 4 You can now mail the return box in a PostNL **post box** ([testalize.me/retourbox](https://testalize.me/retourbox)). After the return box is delivered, you can ask for your result at [testalize.me/uitslag](https://testalize.me/uitslag) within 48 hours, using the code on your code card.